

## AutoFinish Capacity

8/20/24

- **DCPC500:**

- **Zirconia Crown and Bridge:** As many as can fit in a single layer. (Approximately 100)
- **Full Cast Crown and Bridge:** As many as can fit in a single layer. (Approximately 100)
- **Denture:** 6 per barrel.
- **Denture with Metal Exposed:** Only 1 per barrel to avoid damage.
- **Flexible Partial Denture:** As many as can fit in a single layer.
- **Partial Framework:** 6 per barrel.
- **Support Removal:** Model-10 / Denture-6
- **Printed Nightguard, and Splint:** 10 per barrel.
- **Astron Nightguard, and Splint:** 6 per barrel.

- **DCPC250:**

- **Zirconia Crown and Bridge:** As many as can fit in a single layer. (Approximately 50)
- **Full Cast Crown and Bridge:** As many as can fit in a single layer. (Approximately 50)
- **Denture:** 3 per barrel.
- **Denture with Metal Exposed:** Only 1 per barrel to avoid damage.
- **Flexible Partial Denture:** As many as can fit in a single layer.
- **Partial Framework:** 3 per barrel.
- **Printed Nightguard and Splint:** 5 per barrel.
- **Support Removal:** Model-5 / Dentures-3
- **Printed Nightguard, and Splint:** 5 per barrel.
- **Astron Nightguard, and Splint:** 3 per barrel.